

“This is a war. We have won the first battle but we can’t stop there”

Dr. Amy Acton, Director, Ohio Department of Health

These are inspiring words to to make your practice be prepared for the new way of PPE in modern dentistry. **Live Wires Orthodontics** has been on the forefront of 3D printing for over eight years. With this knowledge we are able to provide both the doctor and staff with custom PPE. Life as we know it has definitely changed. With this re-usable PPE it will prove to both your staff and patients that your practice is making their health and well being it's top priority. As we continue battling the corona-19 virus, using custom PPE will help with cross contamination of disposable items, reduce bio hazard waste, and provide reassurance that the doctor and staff will always have adequate PPE available to them. The PPE options are listed as follows:

FULL FACE SHIELDS

\$17.85 ea

- NO quantity Limit
- Available in 8 colors (*clear, red, blue, yellow, orange, green, brown, black*)
- Personal monogram available
- Nylon strap
- Snap closure
- Reusable
- Wipe down sterilization
- Bottle of hand sanitizer included with each full face shield

FULL FACE SHIELDS TO WEAR WITH LOOPS

\$32.85 ea

- No quantity limit
- Available in 8 colors (*clear, red, blue, yellow, orange, green, brown, black*)
- Personal monogram available
- Nylon strap
- Snap closure
- Reusable
- Wipe down sterilization
- Bottle of hand sanitizer included with each full face shield

CUSTOM FACE MASK

\$79.84 ea

- No quantity limit
- Available in 8 colors (*clear, red, blue, yellow, orange, green, brown, black*)
- Monogram mandatory
- Silicone face lining
- Nylon strap
- Snap closure
- Reusable
- Wipe down sterilization
- 1 filter provided
- Bottle of hand sanitizer included with each face mask



FAQs


1. Masks are reusable
2. Masks are not biodegradable
3. Masks can be cleaned with a sterilization product
4. Front face shield can be replaced separately
5. Plastic head band can be replaced separately
6. All materials are FDA approved
7. Masks CANNOT be exposed to extreme heat
8. Masks CANNOT be autoclaved
9. MASKS MUST NOT BE SHARED FROM PERSON TO PERSON
10. PPE equipment cannot be returned once it has been received by purchaser
11. No set up fees
12. No quantity limit
13. Shipping is via the United States Postal Service
14. All materials are manufactured in the United States
15. Masks will take 1-2 weeks, delays may occur with demand



Instructions for scanning a custom face mask. The wearer can be scanned with a standard smart-phone camera.

1. Download BELLUS3D on the iPhone
2. WATCH TUTORIAL
3. Start scanning
4. Export scan as an STL FILE
5. Send STL file to livewirescustommasks@yahoo.com

THE BETTER THE SCANS THE BETTER THE MASK WILL FIT.



LIVE WIRE ORTHODONTICS

904 S. Hazelwood Avenue Youngstown, OH 44509

330.793.9423 | livewiresortho@yahoo.com | livewirescustommasks@yahoo.com